



R&TD Innovative Spontaneous Concepts 2018

# Far-UVC application to spacecraft microbial reduction

Arman Seuylemezian (352N)

Lisa Guan (352N)

David Welch (Columbia University Irving Medical Center)


Manuela Buonnano (Columbia University Irving Medical Center)



**Jet Propulsion Laboratory**  
California Institute of Technology

# Problem Statement

- Current methods for microbial reduction do not mitigate on-going microbial contamination associated with assembly activities, launch




# Project Objective


- Far-UVC is a recently developed technique which utilizes wavelengths at or below 222 nm to inactivate microbes
- The novelty in this wavelength of light is that it's not harmful to humans yet is effective in inactivating microbes.

## Goals:

- Test the efficacy of far-UVC light on inactivating spores of *Bacillus pumilus* SAFR-032 and *Acinetobacter radioresistens* 50v1
- Characterize the effects of varying dose rate on the efficacy of far-UVC
- To investigate far-UVC light as a surface microbial reduction tool for planetary protection purposes




P. Vaishampayan et al., Survival of *Bacillus pumilus* Spores for a Prolonged Period of Time in Real Space Conditions. *Astrobiology* Vol 12, No 5, 2012




K.B. McCoy, I. Derecho, T. Wong, H.M. Tran, T.D. Huynh, M.T. La Duc, K. Venkateswaran, and R. Mogul. *Astrobiology*. Sep 2012. 854-862. <http://doi.org/10.1089/ast.2012.0835>

*Bacillus pumilus* SAFR-032



*Acinetobacter radioresistens* 50v1



# Methods and Approach




Figure 3. Aluminum 6061 coupons were used to seed the bacterial cells/spores into 10 individual spots of 10 ul each containing  $\sim 10^5$  colony forming units (CFU). Spores/cells were recovered by first depositing 200 ul of polyvinyl alcohol and letting it air dry for 1 hour at room temperature then peeling the PVA film and dissolving it into 2 ml of sterile water/PBS respectively.



818-UV/DB low power UV enhanced silicon photodetector with an 843-R optical power meter (Newport, Irvine, Ca)



# Preliminary Data (not from this study)




Figure 1. Measured non-filtered (red) and filtered (blue) UV spectra from our 222 nm KrCl excimer lamp.

Article | [Open Access](#) | Published: 24 June 2020

## Far-UVC light (222 nm) efficiently and safely inactivates airborne human coronaviruses

Manuela Buonanno, David Welch, Igor Shuryak & David J. Brenner

*Scientific Reports* **10**, Article number: 10285 (2020) | [Cite this article](#)

**72k** Accesses | **1815** Altmetric | [Metrics](#)




Figure 2. Bovine corneal tissue exposed to 300 mJ/cm<sup>2</sup> of 222 nm far-UVC light or 254 nm germicidal light. Darkly stained CPD (cyclobutane pyrimidine dimers) positive nuclei are seen in the 254 nm exposed tissue, but the 222 nm exposed tissue shows no increase in CPD over zero-dose controls

# Effects of varying dose rate





Figure 4. Survival of *B. pumilus* SAFR-032 spores exposed to different fluences of 222-nm light delivered at 0.1 or 0.3 mW/cm² (average  $\pm$  SEM).

- Starting population was  $\sim 10^5 - 10^6$  spores per coupon roughly equivalent to M2020 quota
- Demonstrated 3 logs of kill, flattening due to clumping/shielding
- The effects of varying fluence rate are insignificant based on the  $D_{90}$  required to inactivate spores of *B. pumilus* SAFR-032.

# Efficacy of far-UVC in inactivating *Acinetobacter radioresistens* 50v1



- *Acinetobacter* is a non-spore forming microbe, that is hardy and capable of surviving harsh extremes.
- Testing conditions were similar to that of *B. pumilus*, fluence rate was at 0.5 mW/cm²
- There was a significant difference in the dose required to inactivate spores of *B. pumilus* at 0.1 mW/cm² and *A. radioresistens* 50v1 at 0.5 mW/cm² ( $p < 0.005$ ).



# Advantages of far-UVC vs microbial reduction modalities


- Requires minimal intervention by engineers and would therefore have minimal impact on scope of work
- Can be implemented within the assembly environment to continually mitigate microbial bioburden.
- Can be used as a recontamination prevention strategy
  - Launch Pad
  - Payload Fairing during launch
  - Environmental Testing
  - Transport

# Technical Path Forward

1. Link industry standard 254 nm treatment to 222 nm(far-UVC) treatment to ensure acceptance
2. Determine compatibility for spacecraft and facilities materials to establish the potential applications of this technology.
3. Test a wider range of model organisms to demonstrate efficacy across multiple species previously isolated from spacecraft.
4. Establish human- and hardware-safe application protocol
  1. Identify cleanroom friendly lamps
  2. Characterize UV spectrum of candidate lamps
  3. Propose facilities specific implementation plans
5. Perform proof of concept study in an active cleanroom

# Value Added

- Current methods for microbial reduction do not mitigate on-going microbial contamination associated with assembly activities, launch



# Publications



## Abstract

This work aims to investigate far-UVC light at 222 nm as a new microbial reduction tool for planetary protection purposes which could potentially be integrated into the spacecraft assembly process. The major advantage of far-UVC (222 nm) compared to traditional germicidal UVC (254 nm) is the potential for application throughout the spacecraft assembly process in the presence of humans without adverse health effects due to the limited penetration of far-UVC light into biological materials. Testing the efficacy of 222-nm light at inactivating hardy bacterial cells and spores isolated from spacecraft and associated surfaces is a necessary step to evaluate this technology. We assessed survival of *Bacillus pumilus* SAFR-032 and *Acinetobacter radioresistens* 50v1 exposed to 222-nm light on proxy spacecraft surfaces simulated by drying the bacteria on aluminum coupons. The survival fraction of both bacteria followed a single stage decay function up to 60 mJ/cm<sup>2</sup>, revealing similar susceptibility of both species to 222-nm light, which was independent of the exposure rate. Irradiation with far-UVC light at 222 nm is an effective method to decontaminate the proxy spacecraft materials tested in this study.



**Jet Propulsion Laboratory**  
California Institute of Technology

---


[jpl.nasa.gov](https://jpl.nasa.gov)

# Problem Statement

- Our current understanding of the culture based microbial diversity of spacecraft and associated surfaces is based upon 16S rRNA sequences which are limited in resolution.
- The goal of this proposal is to
  - perform taxonomic identification of all the bacterial isolates from the JPL microbial archive collected from spacecraft surfaces from NASA missions and functional genomics using the Whole Genome Sequencing (WGS) data from the most dominant, reoccurring and novel isolates.
- Whole Genome Sequencing will provide
  - More accurate species identifications
  - Functional characterization
  - Subspecies profiling/source tracking
  - A reference database for all future missions utilizing metagenomics



## Genome Encyclopedia of Spacecraft Associated Microbes (GESAM)



## Genome Encyclopedia of Spacecraft Associated Microbes (GESAM)



**Jet Propulsion Laboratory**  
California Institute of Technology






Figure 2. A comprehensive flow chart outlining the GESAM project objectives and tasks.

# What is MALDI-TOF MS

- Matrix Assisted Laser Desorption Ionization Time of Flight Mass Spectrometry
- Generates ions from large molecules while minimizing fragmentation
- Identification based on a protein profile



# Identification Techniques (MALDI-TOF MS)


## Formic Acid Tube Extraction




Inoculate two loops of biomass into 1 ml of 75% ethanol




Vortex for 15 minutes



Pellet isolate and decant ethanol



Add 100 µL 70% formic acid and vortex again for 15 minutes



Add 100 µL 100% acetonitrile and centrifuge x2 minutes



Add 1 µL supernatant to target and overlay with matrix

This document has been reviewed and determined not to contain export controlled technical data.

# MALDI-TOF Background

Goal : Rapid species- level identification of bacterial colonies



Slide provided by: Akemi  
Hinzer

---

# MALDI-TOF Scoring Strategy

Match Type	Score Met	3 Scores	Interpretation
Perfect ID	2.2	✓	Confident in identification
Redo Species	2.2		Identified without confidence
Genus Match	2.0	✓	Confident in genus identification
Redo Genus	2.0		Genus identified without confidence
Taxonomic Group	2.0	✓	Genus identified with mismatched species
Not Reliable	-		All scores were below 2.0
No Peaks			Experiment produced no peaks to match

M2020_2002 (+++)(C)	Bacillus atrophaceus_MER_149B	2.35
M2020_2002 (+++)(C)	Bacillus atrophaceus_MER_149B	2.47
M2020_2002 (+++)(C)	Bacillus atrophaceus_MER_149B	2.53
M2020_2002 (+++)(C)	Bacillus atrophaceus_MER_149B	2.46

M2020_2151 (++)(B)	Bacillus flexus	2.07
M2020_2151 (++)(B)	Bacillus flexus	2.23
M2020_2151 (++)(B)	Bacillus flexus	2.22
M2020_2151 (++)(B)	Bacillus flexus	1.91

M2020_2167 (++)(C)	Bacillus subtilis	2.25
M2020_2167 (-)(C)	not (yet) present	<0
M2020_2167 (++)(C)	Bacillus licheniformis	2.19
M2020_2167 (++)(C)	Bacillus licheniformis	2.21


Slide provided by: Akemi  
Hinzer



# Data Analysis

- Mass spectra of unknown isolates are compared against
  - The Bruker reference database containing:
    - 8540 bacterial strains
  - The JPL in-house curated database containing:
    - 509 bacterial strains

Sample	Reference	Score	Reference	Score
Sample 1	Reference 1	0.95	Reference 2	0.92
Sample 2	Reference 3	0.98	Reference 4	0.95
Sample 3	Reference 5	0.99	Reference 6	0.97
Sample 4	Reference 7	0.96	Reference 8	0.94
Sample 5	Reference 9	0.93	Reference 10	0.91



# 16S rRNA Sanger Sequencing

- Three primer sanger sequencing
  - 27F, 1492R, and 512F
  - More robust, yields higher overlap region than using 2 primers
- Species level identification
  - Although some species are indiscernible using only 16S rRNA
- Still the gold standard for characterizing potentially novel organisms






Figure 1. Three Basic Steps of Automated Sanger Sequencing.

# Genus Level Diversity ~1,200 isolates from spacecraft and associated surfaces



# Species Diversity of 838 strains belonging to the Bacillus Genus



# DNA Extraction, Whole Genome Sequencing and Assembly


- DNA Extraction
  - Manual
    - Filter based extraction kits
  - Automated
    - Magnetic bead based extraction (Maxwell 16 System)
- Quality Check
  - >1 ng/ul is sufficient for sequencing however if extracting from pure cultures >30 ng/ul is reasonable.
  - High molecular weight DNA, running a gel to see the general fragment size is a good idea.



# DNA Extraction, Whole Genome Sequencing and Assembly

## Whole Genome Sequencing

- Illumina paired-end platform (250 bp reads)
- Average of 30-50X coverage (~250 mB per genome raw data)



## Data Analysis

- Quality/length filtering: 0.01% accuracy, >200 bp reads, trim for adapter sequences
- Genome Assembly:
  - <100 contigs
  - L50 (number of contigs needed to reach 50% of the genome): good indication of assembly quality





Welcome to GTDB

# GENOME TAXONOMY DATABASE


317,542 genomes

Release 07-RS207 (8th April 2022)




# How does GTDB-Tk work?

- GTDB-Tk stands for **Genome Taxonomy Database Toolkit** and its workflow consists of 3 steps:
  - Identify
  - Align
  - Classify
- What does this mean for use of the GTDB-Tk Identification database?
  - A more robust method of classification than 16S
  - More parameters used
  - More information including functional genes and applications for Pangenome analysis of isolates




# Taxonomic Identification of 200 whole genome sequenced isolates

## GTDB-Tk Genus Across Missions



Slide provided by Michelle Tran

# GTDB-Tk Species within the Bacillus Genera Across Missions



Slide provided by Michelle Tran

## Next Steps/On-going projects:

- Perform genomic comparisons of strain level variants to identify patterns between species consistently isolated across multiple missions
- Characterize potentially novel organisms
- Identify unique genes present in spacecraft isolated strains in comparison with those from other environments
- Perform metagenomic read mapping against these genomes to identify relative abundance of particular culturable strains across metagenomic datasets.